

COMPLIES WITH REGULATIONS

French Public Health Code article R 1335-1 to R 1335-8

"Any individual or legal entity, be they professionals or private individuals, who produce potentially infectious medical waste, whatever the quantity, is responsible for this waste, from production to disposal."

Local Authorities General Code, articles R2224 -27, circular DGS-VS3/DPPR number 2000/322 of 9/6/00 relating to potentially infectious medical waste being accepted at waste treatment centres.

The local authority is responsible for the disposal of household waste which must be disposed of "in standard conditions without creating a risk to people or the environment."

Decree number 97-1048 of 6th November 1997: Article R44-1a to R44-3, R44-3.

Potentially infectious medical waste must be:

- Separated from other waste from the moment it is produced
- Collected in single-use containers
- Incinerated or pre-processed in accredited disinfection equipment.

7th September 1999 Decree: storage conditions

Specific features of premises:

- They are reserved for the storage of waste and soiled products
 - The stored waste is packed in line with the European Agreement on the International Carriage of Dangerous Goods by Road.
 - The premises are marked, closed, ventilated, lit and equipped with a water inlet and outlet
 - The floor and walls can be washed
- Storage duration between actual production and treatment

- Quantity of more than 100kg per week over 72 hours
- Quantity of less than 100kg per week over 7 days
- Quantity of less than 5kg per month over 3 months

7th September 1999 Decree: monitoring disposal chains

Compulsory documents:

Written agreement between the producer and the waste disposal service provider, a waste tracking slip and a potentially infectious medical waste handover receipt.

PRADAS® Box

Storage centre for potentially infectious and related medical waste



GAP Hygiène Santé

124, rue Saint Sauveur - 06110 LE CANNET
Tél : +33 (0) 493 452 395 • Fax : +33 (0) 493 460 935
Local rate number from within France
0 820 207 240

Visit our website :

www.gap-hygiene-sante.com

STORAGE CENTRE FOR POTENTIALLY INFECTIOUS AND RELATED MEDICAL WASTE



PRADAS® Box

Storage centre for potentially infectious and related medical waste



Ergonomic

Secure

Adaptable

Complies with
regulations

The practical and economical turnkey solution for waste treatment centres, embalming centres, care homes, healthcare centres, ambulance and fire services, blood banks, GDS (Group for Healthy Livestock Breeding) and farms.



PRADAS® Box

Storage centre for potentially infectious and related medical waste

The Pradas® Box is a storage facility which is particularly adapted to storing potentially infectious and related medical waste.

Designed and developed by GAP Hygiène Santé, it complies with regulatory requirements (7th September 1999 decree on storage conditions).



A TURNKEY SOLUTION FOR STORING YOUR POTENTIALLY INFECTIOUS MEDICAL WASTE

Your potentially infectious medical waste is packed in approved containers and stored in the Pradas® Box in optimum health hygiene and safety conditions.

Large storage capacity

TECHNICAL INFORMATION

The Pradas® Box and its technology are protected by French and European patents for which GAP Hygiène Santé holds exclusive rights.

INSTALLATION INFORMATON

- ▶ The Pradas® Box is available in indoor and outdoor versions.
- ▶ **Ready to use with no restrictions:** No need for a construction permit, moveable and adaptable depending on the amount you produce, can be delivered in flat-pack form or assembled.
- ▶ Can be moved and transported on a pallet truck or forklift truck. The cooling unit is easy to detach.

Respect regulations

- Use approved containers which are adapted to the waste in question (needles, blades and sharp waste, liquid waste, soft and solid waste)
- Permanently seal the containers after use
- Store your potentially infectious medical waste in containers or large bulk trolleys in line with the European Agreement on the International Carriage of Dangerous Goods by Road
- You can easily clean your Pradas® Box as often as necessary thanks to the water inlet and outlet system.

STRUCTURE

The insulated chamber is made up of 70mm polyurethane foam sandwich panels with thermal break extrusions:

- ▶ **Inside skin** is made of PVC sheets (optional: stainless steel) with an anti-chafing plinth.
- ▶ **Outside skin** is made of anti-UV lacquered aluminium and is graffiti and salt spray protected.
- ▶ **Reinforced floor** pto support a large bulk trolley (800kg load).
- ▶ **Door** is lockable (refrigerated lorry type) with a door locking system
- ▶ **Feet** 100mm PVC
- ▶ **Electrical power source**
220V single-phase
Outdoor usage version: outside electric cabinet is protected and waterproofed
Indoor usage version: link to mains plug

OPTIONAL FEATURES

- ▶ Automatic waste chain access management and tracking system
- ▶ Door with access code
- ▶ Anti-bacteria system
- ▶ Removable storage shelving and racks – no tools required
- ▶ Reinforced floor for storing large bulk incinerator trolley(s)
- ▶ Detachable aluminium ramp for accessing the large bulk trolleys
- ▶ Full stainless steel interior

OTHER FEATURES

In compliance with the 7th September 1999 decree on storage conditions for potentially infectious medical waste (Article 8) the surface area can be adjusted depending on waste volume and collection frequency.

- ▶ Offers optimum security against the risk of wear and tear and theft:
 - Anti-UV treated sheets are protected against graffiti, aluminium sheets for outside items.
 - Lockable chamber with a refrigerated lorry type lock, cooling unit incorporated inside (combating risk of theft and wear and tear).
- ▶ Specific risks have been identified:
 - Potentially infectious medical waste labelling, yellow colour, Pradas® box logo
 - Fire risk: M2 fire rating
 - Interior locking system: anti-panic lock
 - Certified to comply with EC standards
- ▶ Ventilated, aired, lit – enabling waste to be kept away from heat:
 - Refrigerated (cooling unit) temperature maintained between 8 - 12°C, insulated chamber.
 - Inside waterproof porthole, lighting comes on as door opens
- ▶ Anti-animal intrusion systems: Outside protection grills and anti-intrusion floor drain (rodents etc.)
- ▶ The floor and the walls of the box can be washed
 - Anti-chafing plinths (optional feature: stainless steel interior, anti-bacteria system)
- ▶ Equipped with a water inlet and a washing water outlet
 - Floor drain 150mm x 150mm
 - Water inlet equipped with a HA type shutoff valve and quick connector

OPTION: 30cm stainless steel washing jet tube with 5m hose.

DIMENSIONS

	Pradas® Box for indoor usage	Pradas® Box for outdoor usage
Small volume	1170 x 1220 x 2470 / 215kg (including outside cooling unit)	1220 x 1220 x 2240 / 178kg (inside cooling unit)
Large volume	2130 x 1220 x 2470 / 315kg (including outside cooling unit)	2160 x 1220 x 2240 / 252kg (inside cooling unit)

Other dimensions can be ordered

The Pradas® Box is available for sale or to hire.